



P.O. Box 219 • Batavia, IL 60510  
(800) 452-1261 • Fax (866) 452-1436  
www.flinnsci.com • E-mail: flinn@flinnsci.com  
© 2010 Flinn Scientific, Inc. All Rights Reserved.

**FLINN**  
SCIENTIFIC, INC.  
"Your Safer Source for Science Supplies"

Publication No. 698.00

# Safe Swimming with Sodium

## Introduction

No chemistry class is complete without the spectacular demonstration of alkali metals reacting with water. *Safe Swimming with Sodium* is a novel variation that is much safer to perform than the standard demonstration of simply dropping a small piece of sodium metal into a beaker of water.

## Concepts

- Alkali metals—reaction with water
- Density

## Materials

Sodium metal, Na, 1 small piece	Glass cylinder, approximately 500-mL
Mineral oil, 200 mL	Ring stand and clamp
Phenolphthalein, 1% solution, a few drops	Lithium metal, Li, 1 small piece (optional)
Water, 200 mL	

## Safety Precautions

*Sodium metal is a flammable, corrosive solid and is dangerous when exposed to heat or flame. It will react violently with moist air, water, or any oxidizer. Purchasing pre-cut pieces for performing this demo greatly reduces the potential hazard of the material. Sodium reacts with water to produce flammable hydrogen gas and a solution of sodium hydroxide. Wear chemical splash goggles, chemical-resistant gloves, and a chemical-resistant apron. Please review current Material Safety Data Sheets for additional safety, handling, and disposal information.*

## Procedure

1. Clamp a hydrometer cylinder or large graduated cylinder to a ring stand for support.
2. Add about 200 mL of water to the cylinder followed by a few drops of phenolphthalein solution.
3. Add 200 mL of oil, forming a layer above the water. Tilt the cylinder to reduce mixing at the interface.
4. Drop a piece of sodium, about the size of a kernel of corn, into the cylinder and observe the reaction.

## Disposal

Please consult your current *Flinn Scientific Catalog/Reference Manual* for general guidelines and specific procedures, and review all federal, state and local regulations that may apply, before proceeding. Do not dispose of anything until the sodium has completely reacted. The mineral oil can be stored and reused for future demonstrations and labs. The aqueous solution can be flushed down the drain with excess water according to Flinn Suggested Disposal Method #26a


## Tips

- Mineral oil works best but other hydrocarbon solvents work well. Light mineral oil has low volatility and odor but high viscosity so the sodium will “swim” for a long time. Lighter hydrocarbon mixtures, such as kerosene, turpentine, xylenes, or mineral spirits, have a lower viscosity but a higher volatility and odor, and thus present a greater fire risk if the apparatus is knocked over.

- The colorless water–phenolphthalein layer can be regenerated by the addition of a small amount of dilute acid, such as 1 M HCl. The setup can be used several times during the day.
- Sometimes during the first few reactions, the sodium metal may react very vigorously and briefly melt. If this occurs, the sodium becomes porous and “too light” to sink in the mineral oil. This piece of sodium will no longer swim—try another piece. This sometimes occurs because the mineral oil is wet or becomes wet during the setup.

## Discussion

When added to the cylinder, sodium will sink in the mineral oil until it reaches the interface between the oil and water layers, at which time it reacts with water, forming hydrogen gas and sodium hydroxide, a strong base.



The evolution of hydrogen gas is evident, and hydrogen bubbles adhering to the sodium will carry it into the hydrocarbon layer, temporarily stopping the reaction. The amount of hydrogen and heat evolved is kept under control by this “swimming” behavior, making this demonstration quite safe. The piece of sodium repeatedly dives down to the water–hydrocarbon interface, reacts, then “swims” back up into the hydrocarbon layer until the reaction is complete. During the reaction, the piece of sodium is largely devoid of corrosion, allowing the students to view its gray, metallic appearance. The aqueous layer contains phenolphthalein and turns pink due to the production of a base, sodium hydroxide.

Density is an important physical property that can be used to separate materials or control reactions. Sodium has a density of 0.97 g/mL and sits at the interface of the water and oil layers. Lithium, in contrast, has a density of 0.54 g/cm<sup>3</sup>, and will float on top of the hydrocarbon layer. (Try it!) The interface between two immiscible solvents is an effective site for controlling chemical reactions. Many industrial processes use this concept to react aqueous salts with nonpolar hydrocarbons.

## Connecting to the National Standards

This laboratory activity relates to the following National Science Education Standards (1996):

### *Unifying Concepts and Processes: Grades K–12*

- Systems, order, and organization
- Evidence, models, and explanation

### *Content Standards: Grades 5–8*

- Content Standard B: Physical Science, properties and changes of properties in matter

### *Content Standards: Grades 9–12*

- Content Standard B: Physical Science, structure and properties of matter, chemical reactions

## Acknowledgment

Special thanks to Ken Lyle, St. Johns School, Houston, TX, for bringing this demonstration to our attention.

## Reference

Alexander, M. D., *J. Chem. Ed.* **1992**, *69*, 418.

**Materials for *Safe Swimming with Sodium* are available from Flinn Scientific, Inc.**

Catalog No.	Description
S0329	Sodium, Bottle of 5 small pieces for demonstration
M0064	Mineral Oil, Light, 500 mL
P0115	Phenolphthalein Indicator Solution, 0.5%, 100 mL
AP8599	Hydrometer Cylinder, 600-mL
L0057	Lithium, 2.5 g
AP8916	Safe Swimming with Sodium—Chemical Demonstration Kit

Consult your *Flinn Scientific Catalog/Reference Manual* for current prices.